





**THUASNE®**

## Biflex Self Adjust

|    |  |    |
|----|--|----|
| en | Adjustable compression device for calf and foot .....                    | 6  |
| nl | Verstelbaar compressiehulpmiddel voor kuit en voet .....                 | 6  |
| sv | Justerbar kompressionsanordning för vad och fot .....                    | 7  |
| pl | Regulowana opaska uciskowa na łydkę i stopę .....                        | 8  |
| fr | Dispositif compressif ajustable pour mollet et pied .....                | 9  |
| de | Verstellbare Kompressionsbandage für Wade und Fuß .....                  | 10 |
| it | Dispositivo di compressione regolabile per polpaccio e piede .....       | 11 |
| es | Dispositivo compresivo ajustable para pantorrilla y pie .....            | 12 |
| pt | Dispositivo de compressão ajustável para a barriga da perna e o pé ..... | 13 |
| da | Justerbar kompressionsanordning til læg og fod .....                     | 14 |
| fi | Säädettävä kompressiotuote pohkeeseen ja jalkaan .....                   | 15 |
| no | Justerbar kompresjonsinnretning for legg og fot .....                    | 16 |
| el | Ρυθμιζόμενο σύστημα συμπίεσης κνήμης και άκρου ποδιού .....              | 17 |
| cs | Nastavitelná kompresní pomůcka pro lýtku a chodidlo .....                | 18 |
| lv | Pielāgojams kompresijas izstrādājums apakšstilbam un pēdai .....         | 19 |
| lt | Reguliuojama kompresinė priemonė blauzdai ir pėdai .....                 | 20 |
| et | Reguleeritav kompressiooniseade säärele ja jalale .....                  | 21 |
| sl | Nastavljivi kompresijski pripomoček za meča in stopalo .....             | 22 |
| sk | Nastavitelná kompresívna pomôcka na lýtku a chodidlá .....               | 23 |
| hu | Állítható kompressziós eszköz a lábszárra és a lábfejre .....            | 24 |
| bg | Регулируемо компресиращо изделие за прасец и ходило .....                | 25 |
| ro | Dispozitiv de compresie ajustabil pentru gambă și laba piciorului .....  | 26 |
| hr | Podesiva oprema za kompresiju lista i stopala .....                      | 27 |
| zh | 可调节的小腿或足部压力装置 .....  | 28 |
| ja | ふくらはぎと足用の調節可能な圧迫装置 .....   | 29 |
| ko | 종아리와 발을 위한 조절 가능한 압박 장치 .....  | 30 |
| ar | تجهيرة ضاغطة قابلة للتعديل لريلة الساق والقدم .....                      | 31 |

**Calf - Kuit - Vad - Łydka - Mollet - Wade - Polpaccio - Pantorrilla - Barriga da perna - Læg - Pohje - Legg - Kνήμ - Lýtko - Apakšstilbs - Blauzda - Säär - Meča - Lýtko - Vádlí - Πρασεύ - Gambă - List - 小腿 - ふくらはぎ - 종아리 - ربة الساق**












|   | cB                        | cC                       | cD                         |
|---|---------------------------|--------------------------|----------------------------|
| 1 | 16 - 24 cm<br>6.5 - 9.5"  | 24 - 35 cm<br>9.5 - 14"  | 25 - 36 cm<br>10 - 14"     |
|   | 20 - 28 cm<br>8 - 11"     | 32 - 43 cm<br>12.5 - 17" | 32 - 43 cm<br>12.5 - 17"   |
| 3 | 24 - 32 cm<br>9.5 - 12.5" | 40 - 51 cm<br>15.5 - 20" | 39 - 50 cm<br>15.5 - 19.5" |
|   | 28 - 36 cm<br>11 - 14"    | 48 - 59 cm<br>19 - 23"   | 46 - 57 cm<br>18 - 22.5"   |

| N                             | L                             |
|-------------------------------|-------------------------------|
| $\ell < 36 \text{ cm} / 14''$ | $\ell > 36 \text{ cm} / 14''$ |

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## en **ADJUSTABLE COMPRESSION DEVICE FOR CALF AND FOOT**

### **Description/Destination**

The device comprises a number of adjustable and repositionable overlapping straps which are closed with hook-and-loop fasteners.

Different models are available to fit various parts of the body: foot, calf, knee, thigh, arm and hand.

The calf and thigh models have a tightening system used to control and reproduce the desired level of pressure: low (minimum 20 mmHg), medium (minimum 30 mmHg) and high (minimum 40 mmHg).

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

### **Composition**

#### Components for the calf model:

Textile components: polyamide - elastane - polyurethane - polyester.

Tightening system: polyamide - polyurethane.

#### Components for the foot model:

Polyamide - synthetic rubber (CR) - polyester - elastane.

### **Properties/Mode of action**

The devices provide compression by applying pressure to the limb. They promote venous and lymphatic circulation, help prevent the formation of edema and maintain volume. Inelastic short-stretch material (average pressure between 20 and 50 mmHg).

### **Indications**

Management of edema (venous, lymphatic and caused by injury).

Chronic venous insufficiency (C3 to C6).

Symptomatic treatment of lipedema.

### **Contraindications**

Do not use in case of major dermatological conditions on an affected limb.

Do not use in the event of known allergy to any of the components.

Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6.

Do not use in the event of decompensated heart failure.

Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg).

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of an affected limb.

Do not use in case of extra-anatomical bypass of the affected limb.

Do not use in the event of septic thrombosis.

Do not use on the upper limbs in the following cases:

- Do not use in case of brachial plexus condition.
- Do not use in case of vasculitis of the extremities.

### **Precautions**

Strictly comply with your healthcare professional's prescription and recommendations for use.

If the device seems to be too tight at the recommended tightness, remove it and see a healthcare professional.

The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

For hygiene, security and performance reasons, do not re-use the product for another patient.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the product and seek the advice of a healthcare professional.

Examine the shape of the leg; if there are any deformed or bony areas, use the protective foam bandages or Varico pads to make the shape of the limb uniform.

In the event of a venous ulcer, a Doppler examination is necessary before using the device in order to ensure there is no associated severe arterial disease.

In the event of lipedema, compression can be used to help relieve the symptoms. Seek advice from a healthcare professional for compression and other solutions for managing lipedema.

For some conditions (or situations) such as venous thrombosis, the device should be used in combination with an anticoagulant treatment; consult a healthcare professional for advice.

Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product.

Do not apply the product directly to an open wound without a dressing.

### **Undesirable side-effect**

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness.

Any serious device-related incidents occurring should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

### **Instructions for use/Application**

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

#### Putting the calf model on:

1. Before putting it on, ensure that the hook-and-loop fasteners are rolled up.ⓐ
2. Pull the sleeve ⓑ up to mid-calf height ⓐ
3. If the patient is unable to put the sleeve on, it can be left at the back of the calf.
4. Unroll the bottom strap ⓐ, tighten the hook-and-loop fastener on the other part of the strap (do not fasten to the sleeve) and tighten on the limb according to the desired level of pressure. ⓐ ⓐ
5. Take care not to pull on the adjuster.
6. Repeat with all the other straps from bottom to top (ⓐ, ⓐ and ⓐ) ⓐ
7. Position the top of the product 2 finger breadths below the back of the knee.
8. To remove, undo the straps from top to bottom and roll the hook-and-loop fasteners up again.ⓐ

#### When used with Mobiderm®:

- For the lower limbs, Mobiderm® big blocks pads or bandages should be used.
- Protect the skin with a suitable protective liner (bandage or sock).
- Wrap the Mobiderm® pad around the limb or apply the Mobiderm® bandage in a circular fashion around the limb.
- Position the Biflex Self Adjust device on top.
- Adding a Mobiderm® big blocks pad or bandage increases the circumference by an average of 8 cm. This increase must be taken into account when choosing the size of the device.

#### Using the pressure adjuster ⓐ:

The adjuster comes with 3 levels of pressure: **L** low (light blue line), **M** medium (blue line) and **H** high (dark blue line).

Pressure is reached by pulling on the strap which slides the adjuster until it the outside of the loop is over the desired level of pressure.

When putting it on and wearing it, ensure no large folds appear.

If necessary, the segment pressure can be adjusted by adjusting the loops from bottom to top.

#### Putting the foot model on:

1. Before putting it on, ensure that the hook-and-loop fasteners are rolled up.ⓐ
2. Place the foot in the device with the heel positioned in the hole.ⓐ
3. Unroll the middle strap ⓐ and fasten on the other side while tightening it on the limb.ⓐ
4. Position the bottom strap ⓐ: hold the inner part against the skin, tighten the other part over it and fasten with the hook-and-loop fastener.ⓐ
5. Repeat with the top strap ⓐ.ⓐ
6. To remove, undo the straps from top to bottom and roll the hook-and-loop fasteners up again.ⓐ

#### If the product is too long and extends beyond the toes, it can be trimmed:

1. Determine the appropriate length using the seams on the hook-and-loop fastener ⓐ.
2. Connect the hook-and-loop fastener seam to the corresponding seam at the back and end of the product ⓐ.
3. Cut just above the seams to maintain the shape of the product ⓐ.

#### When used with Mobiderm®:

- For the extremities (hands/feet), Mobiderm® small block pads should be used.
- Protect the skin with a suitable protective liner (bandage or other).
- Cut the Mobiderm® small blocks pad to the desired size and position it over the protective liner at the site of the edema.
- Position the Biflex Self Adjust device on top.

### **Care**

Compression devices are designed for daily use over a 6-month period (i.e. 90 washes). After this time, garments should be replaced to ensure the treatment is effective and technical performance is not affected. Close the hook-and-loop fasteners before washing. Machine washable at 30°C (delicate programme). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

### **Storage**

Store at room temperature, preferably in the original packaging.

### **Disposal**

Dispose of in accordance with local regulations.

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©Thuasne - 2059201 (2025-07)
**THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091 186