**THUASNE**

# Ligastrap® Genu


<b>fr</b>	Genouillère ligamentaire à sanglage fonctionnel.....	6
<b>en</b>	Ligament knee brace with functional strapping.....	7
<b>de</b>	Knieführungsorthese mit funktionellen Zugbändern.....	9
<b>nl</b>	Ligamentaire kniebrace met functionele banden.....	11
<b>it</b>	Ginocchiera con fasciatura funzionale per legamenti.....	12
<b>es</b>	Rodillera ligamentaria con correas funcionales.....	14
<b>pt</b>	Joelheira ligamentar com sistema de correias funcionais.....	16
<b>da</b>	Knæbind til ledbånd med funktionelle stropper.....	17
<b>fi</b>	Polvituki nivelsiteiden hoitoon, toiminnallisuutta edistävät tukisiteet.....	19
<b>sv</b>	Knästöd för ligamentskada med funktionellt remsystem.....	21
<b>el</b>	Επιγονατίδα συνδέσμων με λειτουργικό δέσιμο.....	22
<b>cs</b>	Kolenní vazová ortéza s korekčními pásky.....	24
<b>pl</b>	Orteza stawu kolanowego z funkcjonalnymi pasami dociągowymi.....	26
<b>lv</b>	Ceļa locītavas saišu ortoze ar funkcionālo saiti.....	27
<b>lt</b>	Raiščių antkelis su funkcinė diržų sistema.....	29
<b>et</b>	Funktsionaalse rihmasüsteemiga põlvetugi.....	31
<b>sl</b>	Opornica za koleno s funkcionalnimi trakovi.....	32
<b>sk</b>	Ortéza na podporu kolenných väzov s korekčnými popruhmi.....	34
<b>hu</b>	Ligament térdrögzítő funkcionális pántokkal.....	35
<b>bg</b>	Наколенка за връзките на коляното с функционални ремъци.....	37
<b>ro</b>	Orteză de genunchi ligamentară cu fașe funcționale.....	39
<b>ru</b>	Лигаментарный ортез для коленного сустава с функциональными ремнями.....	41
<b>hr</b>	Ortoza za ligamente koljena s funkcionalnim trakama.....	42
<b>zh</b>	帶有功能性搭扣的韌帶護膝.....	44
<b>ar</b>	داعم للركبة وأربطتها مزود بأحزمة وظيفية.....	47

	1	32 - 34 cm
	2	35 - 37 cm
	3	38 - 41 cm
	4	42 - 45 cm
	5	46 - 49 cm
	6	50 - 53 cm

<b>fr</b>	Circonférence du genou
<b>en</b>	Knee circumference
<b>de</b>	Knieumfang
<b>nl</b>	Omvang van de knie
<b>it</b>	Circonferenza del ginocchio
<b>es</b>	Circunferencia de la rodilla
<b>pt</b>	Circunferência do joelho
<b>da</b>	Omkreds af knæ
<b>fi</b>	Polven ympärys
<b>sv</b>	Omkrets runt knä
<b>el</b>	Περίμετρος του γόνατος
<b>cs</b>	Obvod kolena
<b>pl</b>	Obwód kolana
<b>lv</b>	Celgala apkārtmērs
<b>lt</b>	Kelio apimtis
<b>et</b>	Põlve ümbermõõt
<b>sl</b>	Obseg kolena
<b>sk</b>	Obvod kolena
<b>hu</b>	Térd körfogata
<b>bg</b>	Обиколка на коляното
<b>ro</b>	Circumferința genunchiului
<b>ru</b>	Диаметр окружности колена
<b>hr</b>	Opseg koljena
<b>zh</b>	膝关节周长
<b>ar</b>	محيط الركبة

<b>fr</b>	Stabilisation	Compression	Effet proprioceptif
<b>en</b>	Stabilisation	Compression	Proprioceptive effect
<b>de</b>	Stabilisierung	Kompression	Propriozeptive Wirkung
<b>nl</b>	Stabilisatie	Compressie	Proprioceptieve werking
<b>it</b>	Stabilizzazione	Compressione	Effetto proprioiettivo
<b>es</b>	Estabilización	Compresión	Efecto propioceptivo
<b>pt</b>	Estabilização	Compressão	Efeito propioceptivo
<b>da</b>	Stabilisering	Kompression	Proprioceptisk effekt
<b>fi</b>	Stabiloi	Kompressio	Liikeistinin kohdistuva vaikutus
<b>sv</b>	Stabilisering	Kompression	Proprioceptiv effekt
<b>el</b>	Σταθεροποίηση	Συμπύεση	Αποτέλεσμα ιδιοδεκτικότητας
<b>cs</b>	Stabilizace	Komprese	Proprioceptivní účinek
<b>pl</b>	Stabilizacja	Kompresja	Propriocepcja
<b>lv</b>	Stabilizēšana	Kompresija	Proprioceptīvs efekts
<b>lt</b>	Stabilizavimas	Kompresija	Propriocepinis poveikis
<b>et</b>	Stabiliseerimine	Kompressioon	Proprioseptiivne toime
<b>sl</b>	Stabiliziranje	Kompresija	Proprioceptivni učinek
<b>sk</b>	Stabilizácia	Kompresia	Proprioceptívny účinok
<b>hu</b>	Stabilizálás	Szorítás	Proprioceptív hatás
<b>bg</b>	Стабилизация	Компресия	Проприоцептивен ефект
<b>ro</b>	Stabilizare	Compresie	Efect proprioceptiv
<b>ru</b>	Стабилизация	Компрессия	Проприоцептивный эффект
<b>hr</b>	Stabilizacija	Kompresija	Učinak proprioceptivnog vježbanja
<b>zh</b>	稳定	压紧	本体感
<b>ar</b>	الاستقرار	الضغط	مفعول استقبال الحس العميق

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## en

## LIGAMENT KNEE BRACE WITH FUNCTIONAL STRAPPING

### Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Bilateral product.

Available in 6 sizes and a version made to measure patient.

The device is composed of:

- an anatomically-shaped knitting integrating a light and soft knit at the back of the knee,
- flexible lateral reinforcement,
- a double functional strapping ensuring good lateral support.

### Composition

Textile components: polyamide - polyester - elastodiene - polyurethane - elastane.

Rigid components: polyamide - steel - polyethylene.

### Properties/Mode of action

Stabilisation of the knee ligaments provided by the double strapping system.

An anatomically-shaped knitting for even compression.

## Indications

Joint instability/laxity.  
Knee pain and/or swelling management.  
Sport activities resumption.

## Contraindications

Do not use the product if the diagnosis has not been confirmed.  
Do not use in the event of known allergy to any of the components.  
Do not apply the product in direct contact with broken skin.  
History of venous or lymphatic disorders.  
Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

## Precautions

Verify the product's integrity before every use.  
Do not use the device if it is damaged.  
Choose the appropriate size to fit the patient, referring to the size chart; not applicable to the medical devices made to measure patient, the measures are individually taken.  
It is recommended that a healthcare professional supervises the first application.  
Strictly comply with your healthcare professional's prescription and recommendations for use.  
This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.  
For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.  
In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.  
In the event of any modification in the product's performance, remove it and consult a healthcare professional.  
Do not wear the product in a medical imaging device.  
Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

## Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.  
Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

## Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.  
To avoid damage to the knit fabric, please make sure that the self-gripping pads are always fastened to the designated self-gripping surfaces when putting it on and taking it off.

- Slip the knee brace on the leg with your foot extended, making sure that the patella is properly positioned in the flexible knitted zone provided for this purpose. To make it easier to put on, use the pull tabs at the top of the brace and turn it by approx. 90° before hand.
- To help to put on, a perforation system is available on the strap and on the knee support. The straps Ⓐ have 1 perforation, the straps Ⓑ 2 perforations; be sure to overlap the straps on the corresponding parts on the knee support.
- Pull the straps Ⓐ and fix them at the back on the felt pad at the top of the knee brace (figures 1, 2, 3 and 4).
- Pull the side straps Ⓒ up and fix them half way up on the self adhesive material. Cross them over and place them on the felt pad at the top of the knee brace, in front of where the previous straps were fixed (figures 5 and 6).

• For comfort during wear, make sure the straps are not created.  
If you feel the need to tighten the straps during an activity, make sure you tighten them following the steps indicated in the fitting instructions.

## Care/Maintenance

Close the self-fastening tabs before washing. Machine washable at 40 °C (delicate programme). Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

## Storage

Store at room temperature, preferably in the original packaging.

## Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2023101 (2021-01)