THUASNE®

ROM-R Elbow

fr	Orthèse de coude articulée postopératoire avec contrôle de l'amplitude de mouvement.....	4
en	Hinged post-operative elbow brace with range of motion control.....	4
de	Postoperative Ellenbogenorthese mit Gelenk und Begrenzung des Bewegungsumfangs.....	5
nl	Postoperatieve scharnierende elleboogorthese met regeling van de bewegingsamplitude	5
it	Ortesi di gomito articolata postoperatoria con controllo dell'ampiezza di movimento	6
es	Órtesis de codo articulada postoperatoria con control de la amplitud de movimiento.....	6
pt	Ortótise de cotovelo articulada pós-operatória com controlo da amplitude de movimento	7
da	Artikuleret postoperatorisk ortose til albue med kontrol af bevægelsens udsving	7
fi	Nivelletty kyynärtuki liikelaajuuden hallintaan leikkauksen jälkeen.....	8
sv	Ledad ortos för armbåge för postoperativ användning med kontroll av rörelseomfång	9
el	Μετεγχειρνητική αρθρωτική όρθωση αγκώνα με έλεγχο εύρους κίνησης.....	9
cs	Pooperační kloubová loketní ortéza s kontrolou rozsahu pohybu.....	10
pl	Przegubowa orteza stawu łokciowego z kontrolą zakresu ruchu do stosowania po zabiegach	10
lv	Elkoņa pēcoperācijas ortoze ar šarnīriem, kas ļauj kontrolēt kustības diapazonu	11
lt	Pooperacinis alkūnės įtvaras su jungtimis ir judėjimo amplitudės kontrole	11
et	Operatsioonijärgne liigendatud küünarnukiortoos koos liikumisulatuse kontrolliga	12
sl	Pooperativna zgibna opornica z več možnostmi nadzora gibanja	12
sk	Pooperačná kĺbová ortéza lakťa s kontrolou rozsahu pohybu.....	13
hu	Műtét utáni csuklós könyök ortézis mozgástartomány-szabályozással	13
bg	Постоперативна шарнирна лакътна ортеза с възможност за регулиране на амплитудата на движение.....	14
ro	Orteză postoperatorie articulată pentru cot, cu control al mișcării	14
ru	Послеоперационный шарнирный локтевой ортез с контролем амплитуды движений.....	15
hr	Zglobna ortoza za lakat nakon operacije s kontrolom amplitude pomicanja.....	16
zh	运动范围可控的带铰链术后用肘部支架.....	16
ar	تجهيز مفصلي للمرفق لما بعد الجراحة بميزة مراقبة نطاق الحركة.....	17

Taille Unique / Universal size

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung -
Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso -
Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet -
Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły
użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos -
Kasutusjuhised - Navodila za uporabo - Pokyny na používanie -
Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare -
Инструкция по эксплуатации - Upute za uporabu - 使用指南 -
إرشادات الاستعمال

UK
CA

FDA
REGISTERED

MD

i

i

!

© Thuasne - 6175101 - 2024-09
© Studio Caterin

Importer in EU :
THUASNE
120, rue Marius Auphan
92300 Levallois-Perret
France
www.thuasne.com

Importer in Switzerland :
THUASNE Switzerland SA
Rue du marché 20
1204 Genève
Switzerland

Importer and distributor in the USA :
THUASNE USA / NORTH AMERICA
4615 Shepard Street Bakersfield, CA 93313
Tel. 800.432.3466 - Fax 800.798.2722
www.thuasneusa.com

e-life International Co., Ltd.
7F., No.1, Baosheng Rd., Yonghe Dist.,
New Taipei City 234,
Taiwan (R.O.C.)

EC **REP**

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en

HINGED POST-OPERATIVE ELBOW
BRACE WITH RANGE OF MOTION
CONTROL

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the size chart.

Universal size.

The brace consists of:

- Two lateral adjustable slider bars along the arm and forearm.
- A Range of Motion hinge at the elbow joint that connects the lateral bars.
- 4 adjustable straps to secure it to the patient.
- A strap and handrest (optional).

This brace stabilises, immobilises and helps control elbow joint movement within an appropriate range of motion, helping the patient's healing process.

Composition

Textile components: polyester 25% - polyamide 60% - polypropylene 15%.

Rigid components: aluminum 60% - steel 10% - stainless steel 10% - polyamide 20%.

Properties/Mode of action

The postoperative brace is primarily designed to immobilise and stabilise the elbow joint.

The lateral adjustable slider bars stabilise and immobilise the elbow joint.

The brace can be adjusted in flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) or extension (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°).

Total immobilisation is possible at 0°, 15°, 30°, 45°, 60°, 75° or 90°.

Indications

Postoperative elbow immobilisation or elbow ROM control.

Conservative treatment of elbow dislocation or luxation.

Stable or internally fixed fractures of the distal humerus or proximal radius or ulna.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of unstable or surgical fractures.

Do not use in the event of fractures outside the product's area of action (fracture of the proximal humerus, distal radius or ulna).

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of impaired lymphatic drainage.

Do not use in the event of known allergy to any of the components.

Precautions

Verify the product's integrity before every use.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not use the device if it is damaged.

Do not apply the product in direct contact with broken skin.

Do not use in cases of severe varicose veins or lymphatic drainage disorders.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Do not wear the device in a medical imaging machine.

Do not use in the event of paresthesia.

Do not wear the device for more than 24 hours of continuous use.

It is recommended to wear the device directly on the skin, unless contraindicated.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

Consult a health professional if discomfort or pain occurs.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Additional accessories/spare parts are available to order.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

HEALTHCARE PROFESSIONAL INSTRUCTION
GUIDE:

Please read these instructions with the patient when fitting the brace for the first time and make sure that he/she understands how to position the brace.

PATIENT INSTRUCTIONS GUIDE:

The healthcare professional who initially fitted your brace has made all the necessary adjustments.
The healthcare professional should have also explained how to put on your brace.

INSTRUCTIONS FOR USE

1. Positioning:

- Unfasten the straps and detach the four Quick-clip buckles from the slider bars by rotating the release end upwards until the clip is ejected from the slider bar. ❶
- Place the brace on the lateral side of the arm. ❷
- Ensure that the hinge lock is on the down side.
- Align the hinges with the elbow joint. ❸

2. Length adjustment:

- Pull the pin to slide the upper and lower hinge bars to accommodate arm length and then release the pin. ❹
- Adjust to the same length on the other side.
- To adjust the length, the lower extension (forearm part) must extend to the styloid (wrist).
- The upper extension (upper arm part) should extend to the most comfortable position in the armpit or axilla.
- Starting with the cuffs closest to the elbow, form cuffs around the forearm and biceps.
- Starting at the wrist, click each of the Quick-Clip loops on the sliding upright, making sure that an «audible click» is obtained.
- Tighten the straps for a good fit.
- Pull up to overhang the trapezoid on the opposite side and secure the strap to the appropriate length through the D-ring of the arm slider bar. ❺

3. ROM system adjustment:

To adjust the settings, press the flexion and extension buttons and adjust to the prescribed range of motion. ❻

To lock the range of motion, slide the blue Quick-Lock button into the locked position. ❼

If required, lock ties are supplied with each brace.

Pass a tie through the hole in the extension and flexion buttons to maintain the prescribed range-of-motion settings. ❽

4. Installing the accessories:

Required components: screw (included) - nut (included) - phillips head screwdriver ❶ - flat head screwdriver ❷

- Slide the lower bar down until its hole aligns with the largest hole in the lower bracket. ❸

- Insert the nut into the hole. ②
- Hold the nut in place with one hand and turn it forwards. Use the screw to secure the handrest to the bottom bar. ③
- Use a flat head screwdriver to secure the nut in place. Next, tighten the screw with a Phillips screwdriver to ensure it is secure. ④
- Fasten the second screw as in steps 2 to 4. ⑤
- Fasten the strap as in steps 2 to 4 to complete accessory assembly. ⑥

Care
Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label. Make sure to remove the lateral stays before washing. Hand wash. Do not dry clean. Do not tumble-dry. Wipe pads dry after each use to remove all traces of moisture and allow to air dry. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Dry away from any direct heat source (radiator, sun, etc.). Do not expose the product to extreme temperatures. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage
Store at room temperature, preferably in the original packaging.
Disposal
Dispose of in accordance with local regulations.
Keep these instructions for use.

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MDSS - UK RP Ltd