SpryStep® Pediatric*

fr	Releveur de pied postérieur dynamique.....	3
en	Posterior dynamic ankle-foot orthosis (AFO).....	3
de	Dynamische Peronäusorthese	4
nl	Dynamische enkel-voetorthese posterior.....	5
it	Tutore dinamico posteriore	6
es	Antiequino posterior dinámico	6
pt	Suporte de pé posterior dinâmico	7
da	Dynamisk posterior ankel-fod-skinne	8
fi	Dynaaminen posteriorinen jalannostin.....	9
sv	Dynamisk posterior droppfotsortos	9
el	Δυναμικός οπίσθιος ανυψωτήρας ποδίου.....	10
cs	Dynamická posteriorní ortéza	11
pl	Dynamiczna orteza tylna do podtrzymywania stopy.....	11
lv	Aizmugurējā stiprinājuma dinamiskā krītošās pēdas ortoze.....	12
lt	Dinaminis kulninis pėdos įtvartas.....	13
et	Pöia dünaamiline posterioorne tugi	13
sl	Dinamična opornica za dvig stopala.....	14
sk	Dynamická zadná ortéza na zdvíhanie chodidla.....	15
hu	Hátulról támasztó dinamikus boka-lábfej ortézis.....	15
bg	Динамична ортеза за задна стабилизация на стъпалото	16
ro	Ridicător dinamic posterior pentru picior	17
ru	Динамический ортез-стоподержатель с шиной на задней поверхности голени	18
hr	Ortoza za dinamično podizanje stražnjeg dijela stopala	18
zh	后护壳式动态踝足抬脚器	19
ar	رافع نشيط لمؤخرة القدم.....	20







* Pediatric: Paediatric - Pédiatrique

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso -
 Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání -
 Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie -
 Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 -
 إرشادات الاستعمال

fr	Stabilisation	Correction biomécanique	Restitution d'énergie
en	Stabilisation	Biomechanical correction	Energy return
de	Stabilisierung	Biomechanische Korrektur	Energierückgabe
nl	Stabilisatie	Biomechanische correctie	Energie-teruggave
it	Stabilizzazione	Correzione biomeccanica	Restituzione di energia
es	Estabilización	Corrección biomecánica	Restitución de energía
pt	Estabilização	Correção biomecânica	Restituição de energia
da	Stabilisering	Biomekanisk korrektion	Genoprettelse af energi
fi	Stabiloi	Biomekaaninen korjaus	Energianpalautus
sv	Stabilisering	Biomekanisk korrigering	Frigör energi
el	Σταθεροποίηση	Βιο-μηχανικός έλεγχος	Αποκατάσταση ενέργειας
cs	Stabilizace	Biomechanická korekce	Restituce energie
pl	Stabilizacja	Korekcja biomechaniczna	Odzyskanie energii
lv	Stabilizēšana	Biomehāniskā korekcija	Energijas atgūšana
lt	Stabilizavimas	Mechaninė biologinė korekcija	Energijos atgavimas
et	Stabiliseerimine	Biomehaaniline korrektsioon	Energia taastamine
sl	Stabiliziranje	Popravek biomehanske motnje	Obnovitev energije
sk	Stabilizácia	Biomechanická korekcia	Obnova energie
hu	Stabilizálás	Biomechanikus korrekció	Energia visszanyerése
bg	Стабилизация	Биомеханична корекция	Възвръщане на енергията
ro	Stabilizare	Corecție biomecanică	Recuperarea energiei
ru	Стабилизация	Биомеханическая коррекция	Восстановление энергии
hr	Stabilizacija	Biomehanička korekcija	Obnova energije
zh	稳定	生物力学矫正	体能恢复
ar	الاستقرار	تصحيح ميكانيكي حيوي	استرجاع الطاقة

©

fr

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Trimming value chart				
Footplate length		Product height ❶		
	Original footplate length	Minimum length after trimming*	Original product height	Minimum height after trimming*
XS	14,5 cm	9,5 cm	18 cm	16 cm
	5 ¼"	3 ½"	7"	6 ½"
S	16 cm	11 cm	21 cm	19 cm
	6 ¼"	4 ¼"	8 ¼"	7 ½"
M	17,5 cm	12,5 cm	24 cm	22 cm
	7"	5"	9 ½"	8 ¾"
L	19 cm	14 cm	27 cm	25 cm
	7 ½"	5 ½"	10 ¾"	9 ¾"
XL	20,5 cm	15,5 cm	30 cm	28 cm
	8"	6"	11 ¼"	11"
Calf cuff height = 6 cm (2 ½") (the top and lateral edges of the calf cuff can be trimmed by up to 2 cm (¾")).				
*Do not cut below the minimum value.				

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table. This device is a posterior dynamic ankle-foot orthosis (AFO) that supports and/or stabilizes the foot and ankle while walking.

Composition

Rigid components: glass fibre - aramid fibre.
Textile components: polyamide - elastane.

Properties/Mode of action

The AFO is composed of two parts (rigid and soft part) already assembled. The rigid part, made of composite materials, is positioned under the foot and along the leg to provide stabilisation, biomechanical correction and energy return. The blue translucent areas on the footplate and calf part are trimmable. The textile components (foam pads and strap(s)) protect the leg and allow correct position of the device on the limb.

Indications

- These indications are biomechanical deficits of neurological, traumatic or muscular origin.
- Footdrop.
 - Mild ankle and foot tri-planar instability.
 - Mild knee hyperextension.
 - Hypotonia.
 - High tone.
 - Delayed standing.
- The above may be associated with the following conditions:
- Cerebral Palsy.
 - Poor proprioceptive awareness.
 - Idiopathic toe walking.
 - Spina bifida.
 - Muscular dystrophy.
 - Hypermobility.
 - Selective Dorsal Rhizotomy.
 - Ehlers-Danlos syndrome.

Contraindications

Do not use the product if the diagnosis has not been confirmed.
Do not apply the product in direct contact with broken skin.
Do not use in the event of known allergy to any of the components.
Do not use for patients weighing > 60 kg (130 lbs).
Severe loss of sensation in the lower limb.
Open ulcers of the foot, ankle or lower leg.
Moderate to severe oedema of the affected limb.
Moderate to severe spasticity of the foot and ankle.
Triplanar instability.
Plantarflexion contracture.
Moderate to severe foot deformities.
Moderate to severe ankle instabilities.
Moderate to severe fixed ankle varus or valgus conditions.

Precautions

Verify the product's integrity before every use.
Do not use the product if it is damaged.
The initial fitting and adjustment must be done by a healthcare professional. It is recommended that an adult supervises the application and use of the product by a child.
Strictly comply with your healthcare professional's prescription and recommendations for use.
Regularly check the patient's growth (change in shoe size, floor-to-calf height...).

It is recommended to renew the device (bigger size) if the child's shoe size (foot length) evolves by more than 2 sizes.
Check the condition of the affected limb and the state of the skin daily (with particular attention for patients with sensory deficit).
In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.
For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Do not wear the product in a medical imaging machine.
The systematic use of a sock is recommended when wearing the device.
It is recommended to adequately tighten the device to achieve a good fit on the limb without restricting blood circulation.
Do not expose the product to extreme temperatures.
Avoid excessive pressure on the forefoot area:

- Always place the entire foot on any step or uneven surface. ❹
- Transition sitting/standing position (chair, toilet, car, ...); put the foot flat on the ground before moving to standing. Use any fixed support (armrests, support bar...) to limit the overload of the foot lifter. ❸

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters,etc.) or wounds of various degrees of severity.
Possible risk of venous thrombosis.
Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

The device must be worn in shoes with the following features:
Stiff posterior heel cup sufficiently high to fully encapsulate the foot and the orthosis;

Low heel height;
Shoe closure system: self-fastening tabs or laces.
Trainers or athletic shoes are the best type of shoes to use with the device.
Preparation of the AFO performed by the healthcare professional:
Ensure that the size and the model (right or left) are appropriate for the patient's needs.
Make sure that the patient's footwear complies with the recommendations. The healthcare professional must supervise the fitting of the product and the specific walking conditions of the patient when using the device for the first time.

Sole and calf cuff pattern and trimming: ❸

Note: only the blue translucent areas at the tip of the foot and on the calf part can be cut.

If the shoe has a removable inner sole, use it to draw the cutting outline on the device's sole.

If there is no removable inner sole, position the patient's foot to draw the cutting outline on the device's sole.

Use scissors to cut out the device sole along the traced outline.
The top edge or the side edges of the calf cuff can be trimmed by 2 cm (¾").
Use a buffing cone or fine grade sandpaper to create a smooth finish on trimmed area(s).

Fitting the AFO:

Take the inlay out of the shoe, if removable.
Place the device into the shoe. ❷

Ensure the device's foot plate sits flat on the sole of the shoe and that the posterior heel cup of the shoe is not distorted.

Replace the removable inner sole on top of the device's footplate, unless it makes the shoe too tight. ❶

If there is no removable inner sole then disregard this step.
Loosen the shoe laces and slide the foot inside the shoe. ❶

You may use a shoe horn if required.

Check the fit:

Secure the strap(s):
If a strap is too long, remove the self-fastening tab, trim the strap with scissors and reposition the self-fastening tab.

Tighten the strap(s). ❶

Ensure laces or other tightening mechanisms are firmly fastened. ❶

Ensure comfort of foot and leg, with no impingements, prior to use.
Additional accessories/spare parts are available to order.

Fitting the spare parts:

The spare parts kit contains the following components : foam pad(s), strap(s), self-fastening tab(s).

Remove the textile parts and the self-fastening tab(s) (if damaged) stuck on the rigid part.

Clean the area where the self-fastening tabs were applied.

Replace the self-fastening tab by new one(s) and then position the new foam pad.

If needed, shorten the strap(s): remove the self-fastening tab, trim the strap(s) and replace the tab.

Position the buckle of the strap on the rigid part and follow the fitting instructions.

Care/Maintenance

Product can be washed in accordance with the instructions shown on this leaflet and on the label. If the device comes into contact with water, dry the textile part and wipe the rigid part well with a dry cloth. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Rigid components: wash the rigid part with a moist cloth.

Textile components: the soft part can be fully removed for washing. Replace in the original location before next use. Machine washable at 30 °C (delicate programme). Remove the self-fastening tabs before washing. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not dry clean. Do not tumble-dry. Do not iron. Squeeze out excess water. Dry flat. Dry away from any direct heat source (radiator, sun, etc.).

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

COMMERCIAL WARRANTY AGREEMENT AND WARRANTY LIMITATIONS

Thuasne offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:

- six months for the textile components;
- one year for the rigid components.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:

- any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use,
- damages occurred as part of attempts to modify the product.

Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user or its legal representative (parents, guardian...) to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne entity. Any warranty claim will first be reviewed by Thuasne to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.

To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.

If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product.

It is expressly agreed that this commercial warranty is in addition to the legal warranties binding the entity which sold the product to the user, in accordance with the applicable local legislation in the country of purchase of the product.

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THUASNE DEUTSCHLAND GmbH
THUASNE UK Ltd**