


MAATTABELLEN - CONFECTIE ARTIKELEN

KIT BIFLEX

Zwachtelsysteem voor ulcus cruris

Omvang enkel	Maat	Z-index
<25 cm	1	16811372
25-30 cm	2	16811402
>30 cm	3	16811429




Maat	Hand- omvang (cA)	Pols- omvang (cC)	Elleboog- omvang (cD)	Bovenarm- omvang (cG)
MOBIDERM AUTOFIT ARM				
1	17 - 21 cm	15 - 19 cm	22 - 28 cm	25 - 32 cm
2	17 - 21 cm	16 - 20 cm	24 - 30 cm	29 - 36 cm
3	18 - 22 cm	17 - 21 cm	26 - 32 cm	33 - 40 cm
4	18 - 22 cm	17 - 21 cm	28 - 34 cm	37 - 44 cm
5	19 - 23 cm	18 - 22 cm	30 - 36 cm	41 - 48 cm
6	20 - 24 cm	19 - 23 cm	31 - 37 cm	45 - 52 cm

Universele Z-index: **16187210**

Normaal $\ell = 40 - 45$ cm

Lang $\ell = 45 - 50$ cm




MOBIDERM AUTOFIT ONDERBEEN (AD)

Maat	Enkelomvang (cB)	Kuitomvang (cC)	Onderbeen- omvang (cD)	Kort $\ell D = 38 - 40$ cm	Normaal $\ell D = 40 - 42$ cm	Lang $\ell D = 42 - 44$ cm	Extra lang $\ell D = 44 - 46$ cm
1	16 - 24 cm	24 - 35 cm	25 - 36 cm	16516907	16517091	16517148	16517180
2	20 - 28 cm	32 - 43 cm	32 - 43 cm	16516974	16517105	16517156	16517199
3	24 - 32 cm	40 - 51 cm	39 - 50 cm	16516982	16517113	16517164	16517202
4	28 - 36 cm	48 - 59 cm	46 - 57 cm	16517083	16517121	16517172	16517210

MOBIDERM AUTOFIT BEEN (AG)

Maat	Enkelomvang (cB)	Kuitomvang (cC)	Bovenbeen- omvang (cG)	Kort $\ell G = 67 - 72$ cm	Normaal $\ell G = 72 - 77$ cm	Lang $\ell G = 77 - 82$ cm	Extra Lang $\ell G = 82 - 87$ cm
1	16 - 24 cm	24 - 35 cm	46 - 57 cm	16517229	16517261	16517326	16517369
2	20 - 28 cm	32 - 43 cm	53 - 65 cm	16517237	16517288	16517334	16517377
3	24 - 32 cm	40 - 51 cm	61 - 72 cm	16517245	16517296	16517342	16517385
4	28 - 36 cm	48 - 59 cm	68 - 80 cm	16517253	16517318	16517350	16517393

MOBIDERM STANDARD HANDSCHOEN




Maat	Polsomvang (C)	Handomvang (A)
1	14 - 16 cm	15 - 17 cm
2	16 - 18 cm	17 - 19 cm
3	17 - 19 cm	19 - 21 cm
4	19 - 21 cm	21 - 23 cm
5	20 - 22 cm	23 - 25 cm
6	21 - 23 cm	25 - 27 cm

Z-index:


Zonder vingers: **16187237** Met vingers: **16187229**

FILIFINE MAATTABEL



Maat	Omvang hand (A)
1	13 - 14,5 cm
2	15 - 16,5 cm
3	16,5 - 18 cm
4	18,5 - 20 cm
5	20,5 - 22 cm
6	22,5 - 24 cm

MOBIDERM STANDARD ARM



Maat	Pols- omvang (C)	Elleboog- omvang (E)	Bovenarm omvang (G)
1	14 - 16 cm	22 - 25 cm	25 - 29 cm
1+	14 - 16 cm	24 - 27 cm	28 - 32 cm
2	16 - 18 cm	25 - 28 cm	29 - 33 cm
2+	16 - 18 cm	27 - 30 cm	31 - 35 cm
3	18 - 20 cm	28 - 31 cm	32 - 36 cm
3+	18 - 20 cm	30 - 33 cm	35 - 39 cm
4	20 - 22 cm	31 - 34 cm	35 - 39 cm
4+	20 - 22 cm	33 - 36 cm	38 - 42 cm
5	22 - 24 cm	34 - 37 cm	38 - 42 cm
5+	22 - 24 cm	36 - 39 cm	41 - 45 cm

Universele Z-index: **15778630**

Normaal $\ell = 40 - 45$ cm

Lang $\ell = 45 - 50$ cm

FILIFINE HANDSCHOEN OPEN VINGERS KORT

Maat	Kleur	Z-index (wit)	Z-index (beige)
1	WIT	15639371	15639460
2		15639398	15639479
3		15639401	15639487
4	BEIGE	15639428	15639495
5		15639436	15639509
6		15639444	15639517


MAATTABELLEN - CONFECTIE ARTIKELEN

SKIN TEENKAP



Maat	Omvang (A)	Z-index (zwart)	Z-index (caramel)
1	21 - 22 cm	16307534	16307526
2	23 - 24 cm	16307550	16307607
3	25 - 26 cm	16307577	16307569
4	27 - 28 cm	16307593	16307585

MAATTABEL



Maat	Omvang hand (A)
1	13 - 14,5 cm
2	15 - 16,5 cm
3	17 - 18,5 cm
4	19 - 20,5 cm
5	21 - 22,5 cm
6	23 - 24,5 cm
7	25 - 26,5 cm
8	27 - 28,5 cm

CICATREX HANDSCHOEN LANG GESL. VINGERS



Maat	Z-index (caramel)	Z-index (zwart)
1	15639673	15957411
2	15639681	15957438
3	15639703	15957446
4	15369711	15957454
5	15639738	15957462
6	15639746	15957470
7	15639754	15957489
8	15639762	15957497

CICATREX HANDSCHOEN KORT GESL. VINGERS




Maat	Z-index (caramel)	Z-index (zwart)
1	15639592	15957578
2	15639606	15957586
3	15639614	15957594
4	15639622	15957608
5	15639630	15957616
6	15639649	15957624
7	15639657	15957632
8	15639665	15957640

CICATREX HANDSCHOEN KORT OPEN VINGERS



Maat	Z-index (caramel)	Z-index (zwart)
1	15639770	15639851
2	15639789	15639878
3	15639797	15639886
4	15639800	15639894
5	15639819	15639908
6	15639827	15639916
7	15639835	15639924
8	15639843	15639932

CICATREX BROEK MAATTABEL



Mt	Bekkenomvang (A)	Taille omvang (B)
1	88 - 94 cm	62 - 72 cm
2	95 - 101 cm	67 - 79 cm
3	102 - 108 cm	73 - 85 cm
4	109 - 115 cm	80 - 93 cm
5	116 - 122 cm	90 - 104 cm
6	123 - 129 cm	97 - 113 cm

CICATREX GORDEL



Maat	Z-index
1	15638278
2	15638286
3	15638294
4	15638308
5	15638316
6	15638324

CICATREX GORDEL HOOG



Maat	Z-index
1	15639029
2	15639037
3	15639045
4	15639053
5	15639061
6	15639088

CICATREX BERMUDA



Maat	Z-index (wit)	Z-index (zwart)
1	15637387	15637328
2	15637395	15637336
3	15637409	15637344
4	15637417	15637352
5	15637425	15637360
6	15637433	15637379

CICATREX CAPRI



Maat	Z-index (wit)	Z-index (zwart)
1	15637247	15637328
2	15637255	15637336
3	15637263	15637344
4	15637271	15637352
5	15637298	15637360
6	15637301	15637379


CICATREX LEGGING



Maat	Z-index (wit)	Z-index (zwart)
1	15637514	15637573
2	15637522	15637581
3	15637530	15637603
4	15637549	15637611
5	15637557	15737638
6	15637565	15737646


MAATTABELLEN - CONFECTIE ARTIKELN

CICATREX PARTIEEL HOOFDVERBAND



Maat	Hoofdomvang (B)	Z-index
1	27 - 30 cm	15636860
2	31 - 36 cm	15636879
3	37 - 45 cm	15636887

CICATREX HOOFDVERBAND



Maat	Hoofdomvang (A)	Z-index
0	51 - 52 cm	15636933
1	53 - 55 cm	15636941
2	56 - 57 cm	15636968
3	58 - 59 cm	15636976

CICATREX OORVERBAND




Maat	Hoofdomvang (A)	Z-index
1	< 52 cm	15637220
2	> 52 cm	15637239

CICATREX OPEN HOOFDVERBAND



Maat	Hoofdomvang (A)	Z-index
0	51 - 52 cm	15636895
1	53 - 55 cm	15636909
2	56 - 57 cm	15636917
3	58 - 59 cm	15636925

CICATREX BUIKVERBAND




Maat	Omvang thorax (A)	Z-index
1	70 - 92 cm	15637832
2	93 - 114 cm	15637840

CICATREX VEST ZONDER MOUWEN



Maat	Omvang thorax (A)	Z-index
1	74 - 82 cm	15639096
2	83 - 91 cm	15639118
3	92 - 98 cm	15639126
4	99 - 115 cm	15639134
5	116 - 120 cm	15639142
6	121 - 135 cm	15639150

CICATREX SLEEVE




Maat	Distale omvang (A)	Proximale omvang (B)	Hoogte (AB)	Z-index
1	9 - 11 cm	19 - 23 cm	35 cm	15637859
2	11 - 14 cm	23 - 28 cm		15637867
3	14 - 16 cm	28 - 32 cm	45 cm	15637875
4	16 - 20 cm	32 - 37 cm		15637883
5	20 - 25 cm	37 - 42 cm		15637891
6	25 - 30 cm	42 - 48 cm		15637905
7	30 - 40 cm	48 - 58 cm		15637913
8	40 - 50 cm	58 - 67 cm		15637921

CICATREX SHORT SLEEVE



Maat	Omvang midden van ledemaat (A)	Hoogte	Z-index
1	12 - 14 cm	20 cm	15637948
2	15 - 17 cm		15637956
3	18 - 22 cm	26 cm	15637964
4	23 - 29 cm		15637972
5	30 - 37 cm		15637980
6	38 - 45 cm		15637999
7	46 - 54 cm		15638006
8	55 - 63 cm		15638014

CICATREX LONG SLEEVE



Maat	Distale omvang (A)	Omvang midden van ledemaat (B)	Proximale omvang (C)	Hoogte (AC)	Z-index
1	13 - 15 cm	20 - 25 cm	22 - 28 cm	39 cm	15639525
2	15 - 17 cm	23 - 28 cm	25 - 31 cm		15639533
3	17 - 19 cm	26 - 31 cm	29 - 35 cm		15639541
4	19 - 21 cm	29 - 34 cm	33 - 39 cm		15639568
5	21 - 23 cm	32 - 37 cm	37 - 43 cm		15639576
6	23 - 26 cm	35 - 41 cm	41 - 48 cm		15639584